





Sample Name: CL\_4

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   34
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   84 (F)
Injection Date  : 10/15/2025 12:01:35 PM      Inj       :    7
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.688	BB	0.2145	1.33776e4	984.53802	94.32371
2	11.111	BB	0.1388	491.60617	53.69729	3.46625
3	13.772	BV	0.0936	66.22674	11.09411	0.46696
4	14.134	VB	0.1135	8.44604	1.10094	0.05955
5	16.023	BB	0.0950	180.68610	29.70860	1.27399
6	19.772	BB	0.0801	58.08459	11.28073	0.40955

```

Totals :                      1.41827e4  1091.41969

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_4

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   34
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   84 (F)
Injection Date  : 10/15/2025 12:01:35 PM      Inj       :    7
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*