





Sample Name: CL_ROH_07

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   48
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   87 (F)
Injection Date  : 10/15/2025 8:35:14 PM      Inj       :    3
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.535	BB	0.2749	2.55196	1.28948e-1	0.02134
2	7.682	BB	0.2215	1.17908e4	831.12610	98.61078
3	11.153	BB	0.1597	6.31393	5.77004e-1	0.05281
4	13.821	BV E	0.1053	1.63511	2.23648e-1	0.01368
5	14.126	VB R	0.1091	21.98403	3.01305	0.18386
6	16.031	BB	0.0996	84.20134	12.99481	0.70421
7	19.574	BV	0.0850	34.56100	6.60789	0.28905
8	19.784	VB	0.0894	14.86064	2.49924	0.12428

```

Totals :                      1.19569e4   857.17068

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_ROH_07

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   48
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   87 (F)
Injection Date  : 10/15/2025 8:35:14 PM        Inj       :    3
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***