





Sample Name: CL_Li smore

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   56
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   88 (F)
Injection Date  : 10/16/2025 1:27:09 AM      Inj       :    2
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.532	BB	0.2580	2.24449	1.24684e-1	0.01962
2	7.682	BB	0.2026	1.12587e4	872.08081	98.43218
3	11.139	BB	0.1412	11.06577	1.14055	0.09675
4	13.805	BV	0.0956	3.01665	4.65713e-1	0.02637
5	14.125	VB	0.1033	20.56247	3.02501	0.17977
6	16.030	BB	0.0947	85.92923	14.18592	0.75126
7	19.572	BV	0.0780	37.85122	7.60820	0.33092
8	19.781	VB	0.0855	18.65752	3.32300	0.16312

```

Totals :                      1.14380e4   901.95389
=====

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_Li smore

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   56
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   88 (F)
Injection Date  : 10/16/2025 1:27:09 AM      Inj       :    2
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***