



Sample Name: CL_1

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :    6
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   71 (F)
Injection Date  : 11/5/2025 4:33:24 PM      Inj       :    6
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

| Peak # | RetTime [min] | Type | Width [min] | Area [pA*s] | Height [pA] | Area % |
|--------|---------------|------|-------------|-------------|-------------|----------|
| 1 | 5.616 | BB | 0.2754 | 3.16382 | 1.49198e-1 | 0.00237 |
| 2 | 8.185 | BV R | 0.3343 | 1.31801e5 | 4830.96436 | 98.75634 |
| 3 | 8.913 | VB E | 0.3455 | 24.36190 | 8.62709e-1 | 0.01825 |
| 4 | 13.772 | BV | 0.0770 | 303.24802 | 62.11249 | 0.22722 |
| 5 | 14.155 | VB | 0.0842 | 50.96272 | 9.26511 | 0.03819 |
| 6 | 16.074 | BB | 0.0815 | 1019.23871 | 193.34718 | 0.76370 |
| 7 | 19.788 | BB | 0.0703 | 258.82410 | 55.62361 | 0.19393 |

```

Totals :                      1.33461e5  5152.32465

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_1

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :    6
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    71 (F)
Injection Date  : 11/5/2025 4:33:24 PM        Inj       :    6
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***