





Sample Name: CL\_whiskey\_08

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   41
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   86 (F)
Injection Date  : 10/15/2025 4:18:45 PM        Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

```
=====
                        Area Percent Report
=====
```

```
Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.687	BB	0.2055	1.33165e4	1039.24988	99.07290
2	11.121	BB	0.1394	21.52808	2.33958	0.16017
3	14.123	BB	0.1043	12.34061	1.79238	0.09181
4	16.023	BB	0.0887	83.05923	14.09637	0.61795
5	19.784	BB	0.0828	7.68537	1.42728	0.05718

```
Totals :                      1.34411e4  1058.90548
```

```
=====
                        Summed Peaks Report
=====
```

Signal 1: FID1 A,  
Empty table.

```
=====
                        Final Summed Peaks Report
=====
```

Signal 1: FID1 A,

\*\*\* End of Report \*\*\*