





Sample Name: CL\_ROH\_07

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   53
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   87 (F)
Injection Date  : 10/15/2025 11:37:51 PM      Inj       :    8
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.532	BB	0.2648	2.59646	1.42346e-1	0.02151
2	7.685	BB	0.2137	1.19023e4	880.42804	98.60629
3	11.151	BB	0.1566	6.55641	6.13736e-1	0.05432
4	13.819	BV E	0.1013	1.65668	2.37798e-1	0.01372
5	14.125	VB R	0.1052	22.19137	3.18691	0.18385
6	16.031	BB	0.0968	85.32069	13.67272	0.70685
7	19.574	BV	0.0796	34.92092	6.83608	0.28931
8	19.784	VB	0.0878	14.98630	2.58053	0.12416

```

Totals :                      1.20705e4  907.69816

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_ROH\_07

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   53
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   87 (F)
Injection Date  : 10/15/2025 11:37:51 PM      Inj       :    8
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*