



Sample Name: CL_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   15
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   72 (F)
Injection Date  : 11/5/2025 10:03:31 PM      Inj       :    5
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.616	BV	0.2572	3.16194	1.60966e-1	0.00241
2	6.342	VB	0.3575	3.08904	1.06783e-1	0.00236
3	8.175	BV R	0.3399	1.29098e5	4828.83496	98.55121
4	8.900	VB E	0.0791	54.02125	9.99366	0.04124
5	13.808	BB	0.0826	2.36681	4.41352e-1	0.00181
6	14.148	BB	0.0876	238.33565	43.72483	0.18194
7	16.074	BB	0.0820	1070.35876	201.43661	0.81709
8	19.793	BB	0.0749	526.52588	111.94273	0.40194

```
Totals :                1.30996e5  5196.64190
```

```

=====
                        Summed Peaks Report
=====

```

```
Signal 1: FID1 A,
Empty table.
```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   15
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   72 (F)
Injection Date  : 11/5/2025 10:03:31 PM        Inj       :    5
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***