





Sample Name: CL\_Charlyle

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   64
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   89 (F)
Injection Date  : 10/16/2025 6:19:05 AM      Inj       :    1
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.504	BB	0.3275	1.77298	6.47363e-2	0.01430
2	7.685	BB	0.2150	1.22310e4	897.41022	98.61755
3	11.134	BB	0.1473	16.43850	1.66275	0.13254
4	13.808	BV	0.0970	2.70304	4.09333e-1	0.02179
5	14.131	VB	0.1080	12.93029	1.79562	0.10426
6	16.029	BB	0.0968	120.35027	19.27516	0.97037
7	19.585	BV	0.0807	11.28304	2.16739	0.09097
8	19.795	VB	0.0964	5.97949	9.12649e-1	0.04821

```

Totals :                      1.24025e4    923.69786

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_Charlyle

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   64
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   89 (F)
Injection Date  : 10/16/2025 6:19:05 AM        Inj       :    1
                                           Inj Volume : 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*