



Sample Name: CL_4

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   38
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   74 (F)
Injection Date  : 11/6/2025 12:04:20 PM      Inj       :    8
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.599	BV	0.2930	3.97512	1.72377e-1	0.00302
2	6.335	VB	0.3311	10.01333	4.12646e-1	0.00762
3	8.175	BV	0.3322	1.29900e5	4853.65771	98.80183
4	8.897	VB	0.0691	370.30899	81.30653	0.28166
5	13.770	BV	0.0771	144.68889	29.56468	0.11005
6	14.179	VB	0.0996	5.82854	8.54495e-1	0.00443
7	16.068	BB	0.0846	978.95728	188.47951	0.74459
8	19.785	BB	0.0765	61.52544	12.69631	0.04680

```

Totals :                      1.31476e5  5167.14426

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_4

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   38
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   74 (F)
Injection Date  : 11/6/2025 12:04:20 PM      Inj       :    8
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***