



Sample Name: CL_4

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   39
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   74 (F)
Injection Date  : 11/6/2025 12:41:15 PM      Inj       :    9
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.603	BV	0.2889	3.82510	1.71023e-1	0.00306
2	6.331	VB	0.3273	9.75671	4.13471e-1	0.00782
3	8.149	BV	0.3220	1.23343e5	4765.03662	98.80389
4	8.879	VB	0.0693	351.82498	77.04853	0.28183
5	13.766	BV	0.0764	137.89955	28.51492	0.11046
6	14.172	VB	0.0977	5.40413	8.11676e-1	0.00433
7	16.060	BB	0.0798	927.34460	180.92113	0.74285
8	19.781	BB	0.0755	57.12794	12.00169	0.04576

```
Totals :                1.24837e5  5064.91906
```

```

=====
                        Summed Peaks Report
=====

```

```
Signal 1: FID1 A,
Empty table.
```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_4

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   39
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   74 (F)
Injection Date  : 11/6/2025 12:41:15 PM      Inj       :    9
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***