





Sample Name: CL_Cl arlyl e

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   70
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   89 (F)
Injection Date  : 10/16/2025 9:58:22 AM      Inj       :    7
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whi skey Analysi s FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whi skey Analysi s on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysi s of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capi llary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.687	BB	0.1994	1.21738e4	963.69354	98.63781
2	11.131	BB	0.1318	16.01641	1.80089	0.12977
3	13.806	BV	0.0935	2.67640	4.24761e-1	0.02169
4	14.128	VB	0.0993	12.94352	1.90456	0.10487
5	16.025	BB	0.0918	119.52762	20.56704	0.96846
6	19.582	BV	0.0789	11.22186	2.22361	0.09092
7	19.792	VB	0.0923	5.73527	9.26010e-1	0.04647

```

Totals :                      1.23420e4   991.54042

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_Charlyle

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   70
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   89 (F)
Injection Date  : 10/16/2025 9:58:22 AM      Inj       :    7
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***