



Sample Name: CL\_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   16
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   72 (F)
Injection Date  : 11/5/2025 10:40:06 PM      Inj       :    6
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.611	BV	0.2783	3.54736	1.65345e-1	0.00269
2	6.466	VB	0.3810	3.32182	1.09793e-1	0.00252
3	8.178	BV R	0.3354	1.30155e5	4814.02539	98.55124
4	8.902	VB E	0.0767	52.73133	10.13138	0.03993
5	13.808	BB	0.0823	2.39017	4.47759e-1	0.00181
6	14.149	BB	0.0871	240.39226	44.42070	0.18202
7	16.075	BB	0.0818	1079.53870	203.66779	0.81741
8	19.793	BB	0.0749	531.42847	113.03299	0.40239

```

Totals :                      1.32068e5  5186.00114

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   16
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   72 (F)
Injection Date  : 11/5/2025 10:40:06 PM      Inj       :    6
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

\*\*\* End of Report \*\*\*