



Sample Name: CL_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   11
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   72 (F)
Injection Date  : 11/5/2025 7:37:03 PM      Inj       :    1
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.612	BV	0.2807	3.41401	1.57611e-1	0.00262
2	6.341	VB	0.3623	3.35921	1.11886e-1	0.00257
3	8.174	BV R	0.3296	1.28625e5	4847.21582	98.55215
4	8.899	VB E	0.0837	54.67020	10.01082	0.04189
5	13.809	BV	0.0857	2.53259	4.49497e-1	0.00194
6	14.148	VB	0.0877	239.24881	43.82403	0.18331
7	16.074	BB	0.0817	1063.95447	201.17255	0.81520
8	19.794	BB	0.0747	522.46674	111.36296	0.40031

```

Totals :                      1.30514e5  5214.30517

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   11
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   72 (F)
Injection Date  : 11/5/2025 7:37:03 PM        Inj       :    1
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***