





Sample Name: CL\_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   10
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   82 (F)
Injection Date  : 10/14/2025 9:24:38 PM      Inj       :    1
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

```
=====
                        Area Percent Report
=====
```

```
Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.689	BB	0.2154	1.34129e4	981.48138	96.10108
2	11.122	BB	0.1437	44.46050	4.64540	0.31855
3	13.803	BV	0.0962	2.78890	4.27274e-1	0.01998
4	14.113	VB	0.1063	203.82011	30.38574	1.46033
5	16.033	BB	0.0927	74.32172	11.92036	0.53250
6	19.776	BB	0.0788	218.78474	43.42900	1.56755

```
Totals :                      1.39571e4  1072.28916
=====
```

```
=====
                        Summed Peaks Report
=====
```

Signal 1: FID1 A,  
Empty table.

```
=====
                        Final Summed Peaks Report
=====
```

Signal 1: FID1 A,

Sample Name: CL\_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   10
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   82 (F)
Injection Date  : 10/14/2025 9:24:38 PM      Inj       :    1
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*