



Sample Name: CL_1

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :    1
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    71 (F)
Injection Date  : 11/5/2025 1:29:27 PM        Inj       :    1
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.697	BB	0.2634	2.10666	1.00976e-1	0.00152
2	8.220	BV R	0.3436	1.36611e5	4865.09326	98.77393
3	8.945	VB E	0.0790	4.49901	8.32733e-1	0.00325
4	13.800	BV	0.0818	313.99228	63.32351	0.22703
5	14.184	VB	0.0841	53.01920	9.65669	0.03833
6	16.105	BB	0.0824	1055.08740	197.28873	0.76286
7	19.817	BB	0.0742	267.04208	57.46365	0.19308

```

Totals :                      1.38307e5  5193.75955

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_1

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :    1
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :    71  (F)
Injection Date  : 11/5/2025 1:29:27 PM          Inj       :    1
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***