





```

Acq. Operator   : SYSTEM                               Seq. Line :   50
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                       Location  :   87   (F)
Injection Date  : 10/15/2025  9:48:16 PM              Inj       :    5
                                                    Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                : only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                : vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                : Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

```

Area Percent Report

Sorted By	:	Signal
Mu l t i p l i e r	:	1.0000
D i l u t i o n	:	1.0000
Do not use Multiplier & Dilution Factor with ISTDs		

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.533	BB	0.2478	2.45925	1.38078e-1	0.02057
2	7.685	BB	0.2122	1.17887e4	880.37744	98.60927
3	11.152	BB	0.1488	6.37651	6.15718e-1	0.05334
4	13.820	BV E	0.1000	1.61670	2.35803e-1	0.01352
5	14.125	VB R	0.1051	22.14824	3.18663	0.18526
6	16.031	BB	0.0959	84.23980	13.67090	0.70464
7	19.574	BV	0.0832	34.66222	6.83375	0.28994
8	19.784	VB	0.0869	14.75840	2.57531	0.12345

Totals : 1.19549e4 907.63363

Summed Peaks Report

Signal 1: FID1 A,
Empty table.

Final Summed Peaks Report

Signal 1: FID1 A,

Sample Name: CL_ROH_07

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   50
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location  :   87  (F)
Injection Date  : 10/15/2025 9:48:16 PM        Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***