





Sample Name: CL_ROH_07

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   49
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   87 (F)
Injection Date  : 10/15/2025 9:11:52 PM      Inj       :    4
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.541	BB	0.2713	2.49017	1.30080e-1	0.02033
2	7.682	BB	0.2225	1.20782e4	826.34735	98.60828
3	11.154	BB	0.1598	6.47070	5.72723e-1	0.05283
4	13.823	BV E	0.1056	1.68062	2.28991e-1	0.01372
5	14.128	VB R	0.1163	22.76827	3.00853	0.18588
6	16.034	BB	0.0975	86.44685	13.00716	0.70577
7	19.577	BV	0.0816	35.48829	6.71698	0.28973
8	19.787	VB	0.0944	15.12232	2.50542	0.12346

```

Totals :                      1.22486e4   852.51723

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_ROH_07

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   49
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   87 (F)
Injection Date  : 10/15/2025 9:11:52 PM      Inj       :    4
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***