





Sample Name: CL\_Charlyle

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   66
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   89 (F)
Injection Date  : 10/16/2025 7:32:02 AM      Inj       :    3
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.516	BB	0.3398	1.77102	6.71191e-2	0.01438
2	7.686	BB	0.2014	1.21431e4	948.48462	98.61645
3	11.132	BB	0.1403	16.48627	1.77605	0.13389
4	13.808	BV	0.0936	2.67890	4.24848e-1	0.02176
5	14.130	VB	0.1036	12.85776	1.88478	0.10442
6	16.028	BB	0.0931	119.38286	20.16058	0.96953
7	19.585	BV	0.0792	11.21395	2.21009	0.09107
8	19.795	VB	0.0952	5.97221	9.26951e-1	0.04850

```

Totals :                      1.23134e4   975.93505

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_Charlyle

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   66
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   89 (F)
Injection Date  : 10/16/2025 7:32:02 AM      Inj       :    3
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*