





Sample Name: CL\_whiskey\_08

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   38
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location  :   86 (F)
Injection Date  : 10/15/2025 2:28:23 PM        Inj       :    2
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

```
=====
                        Area Percent Report
=====
```

```
Sorted By          :      Signal
Multiplier         :      1.0000
Dilution           :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.686	BB	0.1956	1.29666e4	1053.59839	99.06714
2	11.120	BB	0.1347	20.93898	2.37966	0.15998
3	14.121	BB	0.1027	12.02438	1.78219	0.09187
4	16.021	BB	0.0878	81.55134	14.04295	0.62307
5	19.784	BB	0.0829	7.58464	1.40694	0.05795

```
Totals :                      1.30887e4  1073.21013
```

```
=====
                        Summed Peaks Report
=====
```

```
Signal 1: FID1 A,
Empty table.
```

```
=====
                        Final Summed Peaks Report
=====
```

Signal 1: FID1 A,

\*\*\* End of Report \*\*\*