



Sample Name: CL_4

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   35
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   74 (F)
Injection Date  : 11/6/2025 10:14:04 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.608	BV	0.2734	3.62683	1.69672e-1	0.00274
2	6.342	VB	0.3365	10.24747	4.14338e-1	0.00775
3	8.183	BV	0.3328	1.30628e5	4810.78027	98.80375
4	8.902	VB	0.0693	370.35135	81.04230	0.28012
5	13.774	BV	0.0777	145.20685	29.37027	0.10983
6	14.183	VB	0.0988	5.83699	8.64456e-1	0.00441
7	16.073	BB	0.0810	984.71228	188.28261	0.74481
8	19.789	BB	0.0766	61.57431	12.69759	0.04657

```

Totals :                      1.32209e5  5123.62151

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_4

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   35
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   74 (F)
Injection Date  : 11/6/2025 10:14:04 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***