





Sample Name: CL\_Cl arlyl e

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   72
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   89 (F)
Injection Date  : 10/16/2025 11:12:17 AM      Inj       :    9
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.676	BB	0.2218	1.21215e4	852.86774	98.63181
2	11.129	BB	0.1477	16.23837	1.58222	0.13213
3	13.805	BV	0.1004	2.64422	3.83832e-1	0.02152
4	14.126	VB	0.1125	12.86309	1.69629	0.10467
5	16.023	BB	0.0986	119.35156	18.66073	0.97116
6	19.581	BV	0.0823	11.17307	2.09156	0.09091
7	19.790	VB	0.0974	5.87507	8.85356e-1	0.04781

```

Totals :                      1.22896e4   878.16773

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_Charlyle

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   72
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   89 (F)
Injection Date  : 10/16/2025 11:12:17 AM      Inj       :    9
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*