





Sample Name: CL\_Li smore

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   59
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   88 (F)
Injection Date  : 10/16/2025 3:16:42 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.539	BB	0.2888	2.31247	1.15886e-1	0.02009
2	7.682	BB	0.2161	1.13285e4	805.19250	98.43277
3	11.144	BB	0.1530	10.92943	1.05343	0.09496
4	13.809	BV	0.0987	3.01800	4.47440e-1	0.02622
5	14.129	VB	0.1087	20.73236	2.85597	0.18014
6	16.033	BB	0.0947	86.54504	13.51112	0.75198
7	19.575	BV	0.0804	38.07029	7.35243	0.33079
8	19.784	VB	0.0872	18.76317	3.25985	0.16303

```

Totals :                      1.15089e4   833.78863

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_Li smore

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   59
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   88 (F)
Injection Date  : 10/16/2025 3:16:42 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*