



Sample Name: CL\_1

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :    7
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   71 (F)
Injection Date  : 11/5/2025 5:10:11 PM      Inj       :    7
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.658	BB	0.2459	2.83059	1.43613e-1	0.00219
2	8.167	BV R	0.3311	1.27740e5	4790.46973	98.75605
3	8.904	VB E	0.3567	22.95637	7.86113e-1	0.01775
4	13.769	BV	0.0768	295.42340	60.72493	0.22839
5	14.150	VB	0.0843	49.46294	8.97262	0.03824
6	16.069	BB	0.0846	987.04437	189.92336	0.76309
7	19.786	BB	0.0745	251.31563	53.75330	0.19429

```

Totals :                      1.29349e5  5104.77366

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_1

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :    7
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    71 (F)
Injection Date  : 11/5/2025 5:10:11 PM        Inj       :    7
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

\*\*\* End of Report \*\*\*