



Sample Name: CL_3

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   25
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   73 (F)
Injection Date  : 11/6/2025 4:09:00 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.608	BB	0.2497	3.14655	1.62644e-1	0.00231
2	8.200	BV	0.3297	1.34338e5	4996.22070	98.67075
3	8.911	VB	0.0711	194.83318	41.25524	0.14310
4	13.783	BB	0.0782	32.65845	6.54306	0.02399
5	14.157	BB	0.0875	535.30591	98.38533	0.39318
6	16.079	BB	0.0819	1039.04773	195.70088	0.76318
7	19.818	BB	0.0950	4.74289	7.37446e-1	0.00348

```

Totals :                      1.36147e5  5339.00530

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_3

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   25
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   73 (F)
Injection Date  : 11/6/2025 4:09:00 AM        Inj       :    5
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***