



Sample Name: CL_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   20
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   72 (F)
Injection Date  : 11/6/2025 1:06:19 AM      Inj       :   10
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.608	BV	0.2535	3.19479	1.62373e-1	0.00243
2	6.345	VB	0.3842	3.24591	1.10017e-1	0.00247
3	8.177	BV R	0.3376	1.29307e5	4810.07275	98.54966
4	8.902	VB E	0.0793	54.61063	10.06719	0.04162
5	13.810	BB	0.0813	2.31261	4.39797e-1	0.00176
6	14.150	BB	0.0872	238.53502	44.06601	0.18180
7	16.075	BB	0.0821	1073.16858	201.55212	0.81790
8	19.794	BB	0.0747	527.93176	112.55257	0.40236

```

Totals :                      1.31210e5  5179.02284

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   20
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   72 (F)
Injection Date  : 11/6/2025 1:06:19 AM        Inj       :   10
                                           Inj Volume: 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***