





Sample Name: CL\_Li smore

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   60
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   88 (F)
Injection Date  : 10/16/2025 3:53:05 AM      Inj       :    6
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.531	BB	0.2633	2.39130	1.27219e-1	0.02101
2	7.683	BB	0.1996	1.12025e4	885.18311	98.43506
3	11.139	BB	0.1378	10.87673	1.15555	0.09557
4	13.806	BV	0.0929	2.94586	4.71211e-1	0.02588
5	14.125	VB	0.1021	20.42307	3.05244	0.17946
6	16.030	BB	0.0935	85.47709	14.34437	0.75108
7	19.573	BV	0.0775	37.64716	7.63659	0.33080
8	19.782	VB	0.0843	18.33855	3.32686	0.16114

```

Totals :                      1.13806e4   915.29736
=====

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_Li smore

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   60
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   88 (F)
Injection Date  : 10/16/2025 3:53:05 AM      Inj       :    6
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*