





Sample Name: CL_4

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   31
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   84 (F)
Injection Date  : 10/15/2025 10:11:39 AM      Inj       :    4
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

| Peak # | RetTime [min] | Type | Width [min] | Area [pA*s] | Height [pA] | Area % |
|--------|---------------|------|-------------|-------------|-------------|----------|
| 1 | 7.693 | BB | 0.2082 | 1.33418e4 | 1022.60986 | 94.32810 |
| 2 | 11.114 | BB | 0.1336 | 489.65396 | 56.24142 | 3.46191 |
| 3 | 13.774 | BV | 0.0916 | 66.15173 | 11.41136 | 0.46770 |
| 4 | 14.136 | VB | 0.1102 | 8.50775 | 1.15092 | 0.06015 |
| 5 | 16.026 | BB | 0.0922 | 180.05327 | 30.79715 | 1.27300 |
| 6 | 19.775 | BB | 0.0792 | 57.87011 | 11.40361 | 0.40915 |

```

Totals :                      1.41441e4  1133.61433

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_4

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   31
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   84 (F)
Injection Date  : 10/15/2025 10:11:39 AM      Inj       :    4
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***