





Sample Name: CL\_1

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :    7
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    81 (F)
Injection Date  : 10/14/2025 7:34:51 PM      Inj       :    7
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.674	BB	0.1963	1.27063e4	1027.45605	96.16187
2	11.117	BB	0.1358	14.42738	1.56098	0.10919
3	13.759	BV	0.0900	155.90887	27.57181	1.17993
4	14.103	VB	0.1000	49.91949	7.66603	0.37779
5	16.017	BB	0.0935	81.83898	13.74500	0.61936
6	19.763	BB	0.0775	205.05447	41.63690	1.55186

```

Totals :                      1.32134e4  1119.63678

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_1

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :    7
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    81 (F)
Injection Date  : 10/14/2025 7:34:51 PM      Inj       :    7
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*