





Sample Name: CL_Li smore

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   55
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   88 (F)
Injection Date  : 10/16/2025 12:50:48 AM      Inj       :    1
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.530	BB	0.2588	2.39028	1.29833e-1	0.02064
2	7.681	BB	0.2028	1.14014e4	882.10504	98.43429
3	11.136	BB	0.1406	11.16256	1.15708	0.09637
4	13.803	BV	0.0933	2.97399	4.73236e-1	0.02568
5	14.122	VB	0.1031	20.79783	3.06567	0.17956
6	16.028	BB	0.0946	86.95331	14.37161	0.75072
7	19.571	BV	0.0779	38.25258	7.70522	0.33026
8	19.780	VB	0.0853	18.82115	3.36495	0.16249

```

Totals :                      1.15827e4   912.37265

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_Li smore

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   55
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location  :   88 (F)
Injection Date  : 10/16/2025 12:50:48 AM      Inj       :    1
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***