





Sample Name: CL\_3

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   26
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   83 (F)
Injection Date  : 10/15/2025 7:09:11 AM        Inj       :    8
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

```
=====
                        Area Percent Report
=====
```

```
Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.692	BB	0.2105	1.35244e4	1021.30902	93.77439
2	11.113	BB	0.1365	214.62044	23.96515	1.48812
3	13.782	BV	0.0898	21.57201	3.60639	0.14957
4	14.114	VB	0.1041	523.31543	80.27568	3.62853
5	16.031	BB	0.0952	77.68636	12.73497	0.53866
6	19.775	BB	0.0793	60.67807	11.93368	0.42073

```
Totals :                      1.44222e4  1153.82490
=====
```

```
=====
                        Summed Peaks Report
=====
```

Signal 1: FID1 A,  
Empty table.

```
=====
                        Final Summed Peaks Report
=====
```

Signal 1: FID1 A,

Sample Name: CL\_3

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   26
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   83 (F)
Injection Date  : 10/15/2025 7:09:11 AM      Inj       :    8
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*