



Sample Name: CL\_3

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   30
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   73 (F)
Injection Date  : 11/6/2025 7:11:25 AM      Inj       :   10
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.610	BB	0.2456	3.13528	1.65097e-1	0.00231
2	8.200	BV	0.3300	1.34167e5	4985.30078	98.67002
3	8.912	VB	0.0709	194.47903	41.35991	0.14303
4	13.784	BV	0.0782	32.74963	6.56648	0.02408
5	14.159	VB	0.0872	534.77472	98.75304	0.39329
6	16.080	BB	0.0850	1038.25220	198.42180	0.76356
7	19.820	BB	0.0968	5.04876	7.67313e-1	0.00371

```
Totals :                1.35975e5  5331.33441
```

```

=====
                        Summed Peaks Report
=====

```

```
Signal 1: FID1 A,
Empty table.
```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_3

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   30
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    73 (F)
Injection Date  : 11/6/2025 7:11:25 AM        Inj       :    10
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

\*\*\* End of Report \*\*\*