





Sample Name: CL_Charlyle

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   65
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   89 (F)
Injection Date  : 10/16/2025 6:55:40 AM      Inj       :    2
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.514	BB	0.3510	1.77316	6.73793e-2	0.01454
2	7.686	BB	0.1981	1.20295e4	960.44293	98.61865
3	11.131	BB	0.1364	16.04339	1.79288	0.13153
4	13.807	BV	0.0920	2.60296	4.21836e-1	0.02134
5	14.129	VB	0.1023	12.65187	1.88499	0.10372
6	16.027	BB	0.0928	118.59728	20.11485	0.97227
7	19.584	BV	0.0791	11.13061	2.19825	0.09125
8	19.794	VB	0.0923	5.69714	9.19727e-1	0.04671

```

Totals :                      1.21980e4    987.84284

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_Charlyle

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   65
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   89 (F)
Injection Date  : 10/16/2025 6:55:40 AM      Inj       :    2
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***