



Sample Name: CL_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   12
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   72 (F)
Injection Date  : 11/5/2025 8:13:42 PM      Inj       :    2
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.610	BV	0.2801	3.48771	1.61394e-1	0.00273
2	6.316	VB	0.3809	3.27766	1.10882e-1	0.00257
3	8.162	BV R	0.3220	1.25871e5	4799.65088	98.55311
4	8.890	VB E	0.0796	53.04143	9.72632	0.04153
5	13.806	BB	0.0829	2.31175	4.28705e-1	0.00181
6	14.142	BB	0.0830	232.14729	43.01427	0.18176
7	16.068	BB	0.0852	1042.05957	198.69977	0.81590
8	19.789	BB	0.0746	511.62399	109.38306	0.40059

```

Totals :                      1.27719e5  5161.17528

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   12
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    72 (F)
Injection Date  : 11/5/2025 8:13:42 PM        Inj       :    2
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***