





Sample Name: CL_4

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   32
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   84 (F)
Injection Date  : 10/15/2025 10:48:14 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.687	BB	0.2012	1.34134e4	1048.74390	94.32328
2	11.108	BB	0.1328	492.85983	57.07157	3.46579
3	13.770	BV	0.0869	66.46574	11.58370	0.46739
4	14.131	VB	0.1091	8.55761	1.17253	0.06018
5	16.022	BB	0.0885	181.18126	30.84526	1.27407
6	19.772	BB	0.0791	58.20569	11.48578	0.40930

```

Totals :                      1.42207e4  1160.90274

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_4

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   32
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   84 (F)
Injection Date  : 10/15/2025 10:48:14 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***