



Sample Name: CL_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   17
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   72 (F)
Injection Date  : 11/5/2025 11:16:39 PM      Inj       :    7
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.608	BV	0.2716	3.36792	1.61312e-1	0.00261
2	6.348	VB	0.3745	3.39376	1.10434e-1	0.00263
3	8.169	BV R	0.3319	1.27176e5	4756.73926	98.55164
4	8.896	VB E	0.0794	53.48690	9.83975	0.04145
5	13.809	BV	0.0830	2.36953	4.38602e-1	0.00184
6	14.147	VB	0.0872	234.67480	43.29849	0.18185
7	16.073	BB	0.0815	1053.96399	200.08955	0.81674
8	19.793	BB	0.0745	517.78217	110.75182	0.40124

```

Totals :                      1.29045e5  5121.42922

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   17
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    72 (F)
Injection Date  : 11/5/2025 11:16:39 PM      Inj       :    7
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***