





Sample Name: CL_ROH_07

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   46
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   87 (F)
Injection Date  : 10/15/2025 7:22:12 PM      Inj       :    1
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.531	BB	0.2757	2.62041	1.31922e-1	0.02184
2	7.676	BB	0.2198	1.18317e4	842.88470	98.60252
3	11.146	BB	0.1560	6.49472	5.92170e-1	0.05413
4	13.817	BV E	0.1057	1.65571	2.25358e-1	0.01380
5	14.121	VB R	0.1097	22.15408	3.01758	0.18463
6	16.027	BB	0.0994	84.51467	13.08463	0.70432
7	18.785	BB	0.1197	7.22409e-1	8.81058e-2	0.00602
8	19.572	BV	0.0803	34.72604	6.72255	0.28940
9	19.782	VB	0.0888	14.80091	2.50947	0.12335

```

Totals :                      1.19994e4    869.25649

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_ROH_07

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   46
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   87 (F)
Injection Date  : 10/15/2025 7:22:12 PM      Inj       :    1
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***