





Sample Name: CL\_whiskey\_08

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   45
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   86 (F)
Injection Date  : 10/15/2025 6:45:46 PM      Inj       :    9
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

```
=====
                        Area Percent Report
=====
```

```
Sorted By          :      Signal
Multiplier         :      1.0000
Dilution           :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.681	BB	0.2036	1.31622e4	1040.34814	99.06614
2	11.118	BB	0.1377	21.16788	2.33544	0.15932
3	14.120	BB	0.1049	12.17101	1.75606	0.09161
4	16.020	BB	0.0895	82.90541	13.91669	0.62399
5	19.784	BB	0.0849	7.83169	1.40888	0.05895

```
Totals :                      1.32863e4  1059.76522
=====
```

```
=====
                        Summed Peaks Report
=====
```

Signal 1: FID1 A,  
Empty table.

```
=====
                        Final Summed Peaks Report
=====
```

Signal 1: FID1 A,

\*\*\* End of Report \*\*\*