





Sample Name: CL_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   13
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   82 (F)
Injection Date  : 10/14/2025 11:14:19 PM      Inj       :    4
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.690	BB	0.2075	1.34869e4	1038.36377	96.10252
2	11.119	BB	0.1370	44.52542	4.94922	0.31727
3	13.802	BV E	0.0965	2.84627	4.34081e-1	0.02028
4	14.112	VB R	0.1002	205.36349	31.46020	1.46334
5	16.031	BB	0.0910	74.54289	12.25652	0.53116
6	19.775	BB	0.0781	219.68961	44.08708	1.56542

```

Totals :                      1.40339e4  1131.55087

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   13
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   82 (F)
Injection Date  : 10/14/2025 11:14:19 PM      Inj       :    4
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***