





Sample Name: CL\_ROH\_07

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   47
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   87 (F)
Injection Date  : 10/15/2025 7:58:50 PM      Inj       :    2
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.534	BB	0.2522	2.52214	1.38657e-1	0.02128
2	7.681	BB	0.2054	1.16885e4	888.70581	98.60934
3	11.148	BB	0.1483	6.36589	6.17328e-1	0.05371
4	13.817	BV E	0.0997	1.59318	2.33286e-1	0.01344
5	14.122	VB R	0.1049	21.91496	3.15907	0.18488
6	16.028	BB	0.0958	83.47682	13.56858	0.70425
7	19.572	BV	0.0785	34.30867	6.84563	0.28944
8	19.783	VB	0.0872	14.65777	2.54629	0.12366

```

Totals :                      1.18533e4   915.81466

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_ROH\_07

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   47
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   87 (F)
Injection Date  : 10/15/2025 7:58:50 PM        Inj       :    2
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*