





Sample Name: CL_ROH_07

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   54
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   87 (F)
Injection Date  : 10/16/2025 12:14:13 AM      Inj       :    9
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.534	BB	0.2560	2.50255	1.35162e-1	0.02129
2	7.683	BB	0.2058	1.15909e4	879.43933	98.61053
3	11.150	BB	0.1486	6.29951	6.09221e-1	0.05359
4	13.819	BV E	0.1000	1.61096	2.34970e-1	0.01371
5	14.124	VB R	0.1041	21.64252	3.15381	0.18413
6	16.030	BB	0.0954	82.62716	13.49302	0.70296
7	19.573	BV	0.0785	34.04220	6.78786	0.28962
8	19.784	VB	0.0872	14.59538	2.53336	0.12417

```

Totals :                      1.17542e4   906.38674
=====

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_ROH_07

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   54
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   87 (F)
Injection Date  : 10/16/2025 12:14:13 AM      Inj       :    9
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***