





Sample Name: CL\_3

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   19
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   83 (F)
Injection Date  : 10/15/2025 2:53:33 AM      Inj       :    1
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.694	BB	0.2047	1.34876e4	1058.07080	93.78772
2	11.113	BB	0.1325	214.08426	24.86351	1.48866
3	13.781	BV	0.0872	21.39234	3.71156	0.14875
4	14.114	VB	0.1011	519.49646	82.86277	3.61237
5	16.030	BB	0.0932	77.47448	13.06117	0.53873
6	19.775	BB	0.0791	60.94194	12.03094	0.42377

```

Totals :                      1.43810e4  1194.60076

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_3

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   19
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   83 (F)
Injection Date  : 10/15/2025 2:53:33 AM        Inj       :    1
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*