





Sample Name: CL_Li smore

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   62
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   88 (F)
Injection Date  : 10/16/2025 5:06:04 AM      Inj       :    8
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.536	BB	0.1845	1.20058	9.43474e-2	0.01039
2	7.684	BB	0.2145	1.13770e4	837.23309	98.44375
3	11.142	BB	0.1461	11.15309	1.10170	0.09651
4	13.808	BV	0.0960	3.03348	4.65883e-1	0.02625
5	14.127	VB	0.1100	20.81106	2.96095	0.18008
6	16.032	BB	0.0923	86.70238	13.99953	0.75023
7	19.574	BV	0.0834	38.14814	7.49724	0.33009
8	19.783	VB	0.0863	18.80468	3.30835	0.16271

```

Totals :                      1.15568e4   866.66110

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_Li smore

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   62
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   88 (F)
Injection Date  : 10/16/2025 5:06:04 AM      Inj       :    8
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***