





Sample Name: CL_3

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   24
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   83 (F)
Injection Date  : 10/15/2025 5:56:09 AM      Inj       :    6
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.683	BB	0.1840	1.29932e4	1148.64807	93.79832
2	11.101	BB	0.1184	205.99211	26.58655	1.48706
3	13.775	BV	0.0873	20.52410	3.78406	0.14816
4	14.106	VB	0.0926	499.50787	84.99709	3.60595
5	16.024	BB	0.0850	74.53918	13.38335	0.53810
6	19.772	BB	0.0773	58.51273	11.91447	0.42240

```

Totals :                      1.38523e4  1289.31359

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_3

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   24
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   83 (F)
Injection Date  : 10/15/2025 5:56:09 AM      Inj       :    6
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***