





Sample Name: CL_1

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :    2
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :    81  (F)
Injection Date  : 10/14/2025 4:31:35 PM      Inj       :    2
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT only
                  .M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe", vol
                  14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                          Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.688	BB	0.2309	1.33499e4	890.65533	96.14385
2	11.138	BB	0.1597	14.96945	1.36758	0.10781
3	13.772	BV	0.0974	164.09955	26.09568	1.18182
4	14.120	VB	0.1119	52.29830	6.93789	0.37664
5	16.033	BB	0.1025	88.61234	13.17680	0.63817
6	19.777	BB	0.0814	215.46185	40.91672	1.55172

```

Totals :                      1.38854e4   979.15000

```

```

=====
                          Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                          Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

*** End of Report ***