





Sample Name: CL_Charlyle

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   67
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   89 (F)
Injection Date  : 10/16/2025 8:08:40 AM      Inj       :    4
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.522	BB	0.3226	1.80357	6.77419e-2	0.01473
2	7.684	BB	0.2005	1.20736e4	948.23273	98.61569
3	11.131	BB	0.1353	16.26211	1.76753	0.13283
4	13.807	BV	0.0940	2.67091	4.21361e-1	0.02182
5	14.129	VB	0.1039	12.73375	1.85910	0.10401
6	16.028	BB	0.0937	118.66006	19.87266	0.96920
7	19.585	BV	0.0799	11.19796	2.18263	0.09146
8	19.795	VB	0.0978	6.15312	9.22300e-1	0.05026

```

Totals :                      1.22430e4   975.32605

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_Charlyle

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   67
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   89 (F)
Injection Date  : 10/16/2025 8:08:40 AM      Inj       :    4
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***