





Sample Name: CL\_Claryllye

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   68
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   89 (F)
Injection Date  : 10/16/2025 8:45:04 AM      Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.686	BB	0.1986	1.20918e4	962.33240	98.62891
2	11.131	BB	0.1380	16.27901	1.79114	0.13278
3	13.806	BV	0.0933	2.69970	4.29494e-1	0.02202
4	14.129	VB	0.1001	12.99465	1.89228	0.10599
5	16.026	BB	0.0919	119.07655	20.46906	0.97127
6	19.584	BV	0.0790	11.15701	2.20680	0.09100
7	19.794	VB	0.0945	5.88646	9.21393e-1	0.04801

```

Totals :                      1.22599e4   990.04256

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_Charlyle

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   68
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   89 (F)
Injection Date  : 10/16/2025 8:45:04 AM        Inj       :    5
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*