



Sample Name: CL_1

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :    4
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   71 (F)
Injection Date  : 11/5/2025 3:19:44 PM        Inj       :    4
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.607	BB	0.2386	2.63650	1.46060e-1	0.00197
2	8.187	BV R	0.3298	1.31997e5	4970.83545	98.75730
3	8.913	VB E	0.3608	25.96371	8.78573e-1	0.01943
4	13.772	BV	0.0771	304.91064	62.25514	0.22813
5	14.154	VB	0.0840	50.94959	9.28775	0.03812
6	16.074	BB	0.0813	1017.51746	193.69237	0.76129
7	19.789	BB	0.0745	258.99100	55.49227	0.19377

```

Totals :                      1.33658e5  5292.58761
=====

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_1

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :    4
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   71  (F)
Injection Date  : 11/5/2025 3:19:44 PM         Inj       :    4
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***