



Sample Name: CL_3

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   23
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   73 (F)
Injection Date  : 11/6/2025 2:55:55 AM        Inj       :    3
                                           Inj Volume: 0.5 µl
Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Signal
Multiplier          :      1.0000
Dilution            :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	5.609	BB	0.2286	2.79023	1.59113e-1	0.00209
2	8.187	BV	0.3331	1.31657e5	4906.28662	98.66995
3	8.904	VB	0.0716	191.50667	40.21383	0.14352
4	13.781	BB	0.0780	32.07833	6.45570	0.02404
5	14.153	BB	0.0827	524.25909	97.62039	0.39291
6	16.074	BB	0.0811	1019.05756	194.71526	0.76373
7	19.816	BB	0.0972	5.00883	7.19802e-1	0.00375

```

Totals :                      1.33431e5  5246.17071

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_3

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   23
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :    73 (F)
Injection Date  : 11/6/2025 2:55:55 AM         Inj       :    3
                                           Inj Volume : 0.5 µl

Sequence File   : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\FA25 chem.S
Method          : D:\105 F25\Data\chem\2025-11-05 13-27-22 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".

=====
```

*** End of Report ***