





Sample Name: CL_Evan

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   73
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   90 (F)
Injection Date  : 10/16/2025 11:49:26 AM      Inj       :    1
                                           Inj Volume: 0.5 µl

Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	6.514	BB	0.2926	2.84730	1.31484e-1	0.02291
2	7.677	BB	0.2045	1.22531e4	962.79486	98.60243
3	11.153	BB	0.1501	3.58444	3.42552e-1	0.02884
4	13.782	BV	0.0895	7.85169	1.31801	0.06318
5	14.117	VB	0.1036	15.70869	2.30137	0.12641
6	16.019	BB	0.0896	87.46686	14.65515	0.70386
7	19.562	BV	0.0782	42.51168	8.52149	0.34210
8	19.774	VB	0.0877	13.70310	2.35984	0.11027

```

Totals :                      1.24268e4   992.42476

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL_Evan

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   73
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location :   90 (F)
Injection Date  : 10/16/2025 11:49:26 AM      Inj       :    1
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

*** End of Report ***