





Sample Name: CL\_2

```

=====
Acq. Operator   : SYSTEM                      Seq. Line :   12
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown              Location  :   82 (F)
Injection Date  : 10/14/2025 10:37:46 PM      Inj       :    3
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====

```

```

=====
                        Area Percent Report
=====

```

```

Sorted By      :      Signal
Multiplier     :      1.0000
Dilution      :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Type	Width [min]	Area [pA*s]	Height [pA]	Area %
1	7.694	BB	0.1993	1.33523e4	1057.44409	96.10786
2	11.120	BB	0.1347	44.12978	5.01428	0.31764
3	13.801	BV	0.0921	2.69926	4.36595e-1	0.01943
4	14.112	VB	0.0984	202.61313	31.75763	1.45837
5	16.031	BB	0.0941	74.02426	12.32097	0.53281
6	19.774	BB	0.0777	217.27188	43.93592	1.56389

```

Totals :                      1.38931e4  1150.90948

```

```

=====
                        Summed Peaks Report
=====

```

```

Signal 1: FID1 A,
Empty table.

```

```

=====
                        Final Summed Peaks Report
=====

```

Signal 1: FID1 A,

Sample Name: CL\_2

```
=====
Acq. Operator   : SYSTEM                      Seq. Line :   12
Sample Operator : SYSTEM
Acq. Instrument : Charlie Brown                Location :   82 (F)
Injection Date  : 10/14/2025 10:37:46 PM      Inj       :    3
                                           Inj Volume: 0.5 µl
Method          : D:\105 F25\Data\chem\2025-10-14 15-51-51 FA25 chem\Whiskey Analysis FRONT
                  only.M (Sequence Method)
Last changed    : 9/4/2025 3:10:57 PM by SYSTEM
Method Info     : Whiskey Analysis on DUAL SBP-20 columns. Based on "The Reporter Europe",
                  vol 14, November 2004 International Issue, p Supelco doc."The Analysis of
                  Alcoholic Beverages on the 30m x 0.25mm ID, 1.0um SPB-20 Capillary Column".
=====
```

\*\*\* End of Report \*\*\*